
CITIZEN

INSTRUCTION MANUAL FOR DIGITAL EAR THERMOMETER

Model: **CTD504**



(Antibacterial Model)

English
Español

**Instruction Manual for
Digital Ear Thermometer
Model: CTD504 (Antibacterial Model)**

PARTS NAMES / Display Symbols

Accessory:
One +3V CR2032 Lithium Battery

Note: Any artificial changes to the original goods may result in measurement error.

Symbol	Description
①	Measuring Result. Displays the result in Fahrenheit or Celsius.
②	Low Battery Symbol. Flashes to alert the user to change the battery.
③	Too High Ambient Temperature. Refer to “Troubleshooting” for details.
④	Too Low Ambient Temperature. Refer to “Troubleshooting” for details.
⑤	Ready to Use Symbol. The CTD504 is ready to use.
⑥	Commutable Temperature Unit. Our default unit is Celsius. By adjust settings, user can have Celsius or Fahrenheit unit as they see fit. See how to switch from $^{\circ}\text{C}$ to $^{\circ}\text{F}$ (from $^{\circ}\text{C}$ to $^{\circ}\text{F}$) page for details.
⑦	Memory Mode Symbol. With this symbol, the LCD screen displays the previous result.

PRECAUTIONS (Intended Use)

1. The CTD504 is a electronic thermometer that uses an infrared sensor to detect the body temperature of people of any age from the auditory canal in a home setting.
2. Do not use the CTD504 as a replacement for medical advice. Consult your doctor if you have any concerns.
3. Do not immerse CTD504 in water, or any other liquid.
4. Do not drop the CTD504, and avoid any strong impact.
5. Do not use the CTD504 on any body part other than your ear.
6. Do not move or shake your body while measuring.
7. Do not use for any purpose other than the intended use. The manufacturer is not liable for misuse.
8. The CTD504 can measure up to six times a day. To know your normal body temperature range, we suggest you measure regularly.
9. If any damage occurs to the sensor inside the probe tip or to the unit itself, return the product to your local distributor.
10. Keep the CTD504 away from children. Do not let them use the CTD504 on their own. They may choke on the battery or tiny component parts by accident.
11. Install a new CR2032 lithium battery when the power is low. Use of any other battery may lead to inaccurate results or malfunction.

12. Without regular recalibration, the results may not be accurate. Every CTD504 is calibrated right after production. See “Warranty and Recalibration” for details.
13. Keep the CTD504 away from sunlight and dusty places.

Important Notice:

The CTD504 does not require a probe cover during use. It is important to clean and disinfect the CTD504 before and after each use. See “Cleaning and Disinfection” for details.

For measurement:

- When asleep, pressure on the ear may cause higher-than-normal temperature. Wait a few minutes, and then measure again.
- Clean your external ear before measuring. To ensure an accurate result, keep your ear free from obstructions or excess earwax.
- The measuring results from the left and right ears may differ slightly. We recommend measuring in the same ear.
- Before measuring, the user and the CTD504 should be in a steady room temperature for at least 30 minutes. The user should also rest for 30 minutes after exercising.
- Do not modify this equipment without authorization of the manufacturer.
- The patient is an intended operator.
- CTD504 is an adjusted mode (ear mode) clinical thermometer.

HOW TO OPERATE

1. Follow steps in “Cleaning and Disinfection” section to clean CTD504 up.
2. Press the  button. The LCD screen displays all digital segments. 
3. If the CTD504 has saved past results, the LCD screen automatically displays the latest result. If not, skip to next step.
4. When  or  flashes with a short beep, the CTD504 is ready to use.
 1. To improve accuracy, follow the illustrations below. For adults, pull the ear up and back. For infants and children, pull the ear back.

2. Fit the probe tip into the ear canal, and then press the **START** button. Release the button and wait a few seconds until you hear a long beep. Each measurement takes about 1 ~ 3 seconds. 
3. Take the CTD504 out of the ear canal. The LCD displays the latest measuring result. Without any operation, the CTD504 automatically shuts off after 1 minute.

Note: **Before re-measurement, wait 1 – 3 minutes.**

HOW TO SWITCH FROM °C TO °F (FROM °F TO °C).

1. In power-off mode, press and hold the **START** button for 6 seconds. When °C starts flashing, release the button.
NOTE: Do not attempt to turn CTD504 off in the process.
2. Press the **START** button again to switch from °C to °F.
3. Wait until the CTD504 automatically shuts off and restarts. It is now ready to use.

Important Notice before Measuring

- In sleep, pressure on the ear may cause higher-than-normal temperature. Wait several minutes, and then measure again.
- Make your external ear clean before measuring. To secure an accurate result, keep your ear free from obstructions or excess earwax.
- The measuring results took from left or right ears differ slightly. We suggest you measure in the same ear.
- Before measuring, users and CTD504 should be in a steady room for at least 30 minutes. User should also take a 30-minutes-rest after exercising beforehand.

MEMORY STORAGE

CTD504 can store up to 12 measuring results at once. When full, CTD504 overrides the oldest data with new one. Next time when you start CTD504 up, LCD screen automatically displays last measuring result, with flickering **M** Symbol.

RECALL MEMORY

1. In power-off mode, press and release the **START** button to enter Memory Recall Mode. Press the **START** button again. The LCD screen displays the first memory set number with the latest measuring result.
2. Press the **START** button repeatedly to scroll through all 12 stored results in a loop.

ERASE MEMORY

1. While holding the **START** button, press and release the **I/O** button. The CTD504 starts up.
2. When the LCD displays the previous measuring result release the **START** button.
3. To check memory erasing, turn the CTD504 off, and enter Memory Recall Mode again.
Once deleted, your data can not be restored.


BATTERY REPLACEMENT

- flashes: You can still measure, but the accuracy may not be impaired.
- lights or no reaction at all: Change the battery as soon as possible.

Note: Change the battery away from small children and any source of heat.

Small children may choke on battery, or tiny component parts by accident.

1. Press the upper part of the battery cap to unlock it.
2. Insert the tip of a tool into the gap at the side of the battery compartment. Remove the worn-out battery.
3. Insert a new CR2032 lithium battery into the CTD504, according to the polarities in the figure. Push it firmly into place with a "click", then replace the cap.

NOTE: Please dispose of worn-out device and batteries per local regulations.

TIPS FOR NEW USERS

We already put one CR2032 lithium battery inside CTD504. Open battery cap as the figure.


CLEANING INSTRUCTIONS

The CTD504 does not require a probe cover during use. Clean and disinfect the CTD504 before and after each use. If you wish to share the CTD504 with others, follow the instructions below.

◆ Preparation

1. Cotton bud: Available in supermarkets.
2. 70 % concentration alcohol swab: Available in supermarkets.

Note: Dispose of used cotton buds and alcohol swabs in accordance with local regulations.

◆ Cleaning : Probe area

Use a clean cotton bud to wipe off any visible dirt or ear wax. Gently wipe the area clean.


Note: The most delicate and sensitive part of the CTD504 is the sensor inside the probe tip. Poking it with force may cause damage to the CTD504. When cleaning, take extra care.

Hold the CTD504 facing down to prevent water getting into the CTD504. Use an alcohol swab to clean the probe area in a circular motion for 15 seconds (around 20 times).

Afterwards, use the same alcohol swab to clean exterior area back and forth for another 15 seconds (around 20 times).

DISINFECTION (Probe & Exterior Area)

1. Hold the CTD504 facing down to prevent water getting into the CTD504. Use an alcohol swab to clean the probe area in a circular motion for 15 seconds (around 20 times).
 2. Then, use the same alcohol swab to wipe the exterior area back and forth for another 15 seconds (around 20 times).
 3. Wait at least 30 seconds until the CTD504 is completely dry.
- Continue to set units on the next page, and start measuring. Repeat the steps above when measurement is done.

NOTE:

Do not use detergent or any strong chemicals on the CTD504.

Be sure to clean CTD504 up before disinfection. Please wash your hands with soap and water before disinfection.

◆ Signs of Deterioration

During Cleaning and Disinfection, beware of following signs of deterioration:

- Glossy surface of CTD504 blurs.
- No reaction at all toward operating.
- Indelible water stain on the LCD screen.
- Unable to measure.

If any signs listed above occur, please contact your local distributor.

FEVER ALERT

The CTD504 has a built-in "Fever Alert" feature.

If you have a fever between 37.5 °C - 43.0 °C (99.5 °F - 109.4 °F), the "Fever Alert" function will inform you by a flashing result, and a series of beeps.

NOTE:

Do not use the CTD504 as a replacement for medical advice. Consult your doctor if you have any concerns.

MAINTENANCE and STORAGE

◆ Storage

1. If you will not use the CTD504 for over 2 months, remove the battery from the compartment.
2. Always store the CTD504 in its container after use.
3. Keep the CTD504 away from:
 - Extreme temperatures
 - Overly humid or dry environment
 - Direct sunlight
 - Dusty places.
4. For details of environmental conditions, refer to “Specifications”. Using the CTD504 outside the stated range may lead to inaccurate results or malfunction.

◆ Maintenance

1. Do not disassemble or change any parts of the CTD504. Such action may reduce its accuracy.
2. Do not drop the CTD504, and avoid any strong impact.

ABOUT THE EAR THERMOMETER


Normal body temperature differs by body parts, time and age. Here are two tables for your reference:

1. Normal Body Temperature Range

Body Site	Normal Range
Oral	35.5 ~ 37.5 °C 95.9 ~ 99.5 °F
Axillary	34.7 ~ 37.3 °C 94.5 ~ 99.1 °F
Rectal	36.8 ~ 38.0 °C 98.2 ~ 100.4 °F
Ear	35.8 ~ 37.7 °C 96.4 ~ 99.9 °F

Chamberlain, J.M., Terndrup, T.E., The ears have it. Using and interpreting tympanic thermometers., Advance PA, February 1995.

2. Figure: Diurnal Variations in Body Temperature

Source: 1911, Encyclopedia Britannica, 11th edition, Volume 2, Part 1, Slice 1, p95.

TROUBLESHOOTING

LCD screen	Description & Solution
	<p> flashes</p> <p>Battery power is low, but the CTD504 can still measure.</p> <p>→ Remove the worn-out battery. Insert a new CR2032 battery as soon as possible. Without enough energy, the CTD504 cannot function as it should.</p>
	<p> lights</p> <p>Battery power is too low to measure.</p> <p>→ Remove the worn-out battery. Insert a new CR2032 battery as soon as possible.</p>
	<p>Too high ambient temperature</p> <p>Ambient temperature is higher than 40.0 °C (104.0 °F)</p> <p>→ Leave the CTD504 in a room around 10.0 °C – 40.0 °C (50.0 °F – 104.0 °F) for 30 minutes. Then measure again.</p>
	<p>Too low ambient temperature</p> <p>Ambient temperature is lower than 10.0 °C (50.0 °F)</p> <p>→ Leave the CTD504 in a room around 10.0 °C – 40.0 °C (50.0 °F – 104.0 °F) for 30 minutes. Then measure again.</p>

LCD screen	Description & Solution
	<p>Outside displayed temperature range I</p> <p>Result is higher than 43.0 °C (109.4 °F).</p> <p>→ The user should rest in a steady room temperature for at least 30 minutes. Then measure again.</p>
	<p>Outside displayed temperature range II</p> <p>Result is lower than 34.0 °C (93.2 °F).</p> <p>→ The user should rest in a steady room temperature for at least 30 minutes. Then measure again.</p>
	<p>Error symbol I</p> <p>Disturbance occurred while measuring.</p> <p>→ Repeat the measuring steps again. If the error persists, return the CTD504 to your local distributor.</p>
	<p>Error symbol II</p> <p>Memory error occurred while measuring.</p> <p>→ Repeat the measuring steps again. If the error persists, return the CTD504 to your local distributor.</p>

OTHER INFORMATION

This product conforms to the provisions of the EC directive 93/42/EEC (Medical Device Directive). And also complies with following standards (included but not limited):

Safety standard:

EN 60601-1 Medical electrical equipment Part 1: General requirements for basic safety and essential performance

EMC standard:

EN 60601-1-2 Medical electrical equipment Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests

Performance standards:

EN 12470-5 Clinical thermometers - Part 5: Performance of infra-red ear thermometers (with maximum device)

Follow instructions for use.

Temperature limitation

IP22

Degrees of protection provided by enclosures

**TYPE BF APPLIED PART**

- Internally powered equipment
- Not suitable for use in presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide
- Continuous operation with short-time loading

To avoid inaccurate results caused by electromagnetic interference between electrical and electronic equipments, do not use the device near a mobile phone or microwave oven. At least keep a maximum output power of 2 w yields and a distance 3.3m away from this equipment.

Discard the used product to the recycling collection point according to local regulations.

Manufacturer:

HEALTH & LIFE Co., Ltd.

9F, No. 186, Jian Yi Road, Zhonghe District, New Taipei City, Taiwan

www.healthandlife.com.tw

EC REP Authorized Representative in the European Community:

EMERGO EUROPE

Prinsessegracht 20, 2514 AP The Hague, The Netherlands

Date of Manufacture

Factory:

Health & Life (Suzhou) Co., Ltd.

SN Serial Number

No.1428 Xiang Jiang Road, Suzhou New District,

Suzhou City 215129, Jiangsu Province, China

LOT Batch code

Warranty & Recalibration

Warranty For One Year from Manufacturing Date

CITIZEN SYSTEMS JAPAN Co., Ltd. grants a 1 year warranty on the product commencing on the date of purchase. Within the limited time, we will eliminate any defects in the appliance resulting from faults in the material or workmanship free of charge.

Please note that this service does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorized accessory; alteration to the product; improper installation; unauthorized repairs or modification; improper use of electrical/power supply; loss of power; dropped product; malfunction or damage of an operation part due to failure to provide the manufacturer's recommended maintenance; transportation damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of the importers or distributors.

Recalibration Notice

Every CTD504 is properly calibrated right after production.

We recommend that USERS bring the CTD504 to certain laboratories every 2 years for proper recalibration.

The CTD504 does not need periodic recalibration, if the user handles it according to this manual.

SPECIFICATIONS

Measurement method	Infrared
Dimension (L x W x H)	120 x 39 x 52 mm / 4.73 x 1.54 x 2.03 inch
Weight (without battery)	66.8 g / 2.36 oz
Measurement range	34.0 °C – 43.0 °C (93.2 °F – 109.4 °F)
Low temperature display	Temperature < 34.0 °C (93.2 °F): display: L °C (L °F)
High temperature display	Temperature > 43.0 °C (109.4 °F): display: H °C (H °F)
Display resolution	0.1 °C or 0.1 °F
Maximum Permissible Range	35.5 °C ~ 42.0 °C (95.9 °F ~ 107.6 °F) : ± 0.2 °C / ± 0.4 °F Other Range : ± 0.3 °C / ± 0.5 °F
Display	Liquid Crystal Display
Clinical accuracy (repeatability)	± 0.3 °C / ± 0.5 °F
Memory sets	12
Auto power off	1 minute
Battery	CR2032 lithium battery (3V) x 1
Battery life	1500 times
Operating ambient conditions	Temperature: 10.0 °C ~ 40.0 °C (50.0 °F ~ 104.0 °F) Humidity: 15 % ~ 95 % R.H. Atmospheric pressure: 700 hPa ~ 1060 hPa
Storage / Transportation Conditions	Temperature: -25 °C ~ 55 °C (-13.0 °F ~ 131.0 °F) Humidity: ≤ 95 % R.H.
Safety classification	 Type BF applied part
Product Life	5 Years
Accessory	CR2032 lithium battery (3V) x 1

* The contents and specifications in this manual are subject to change for improvement without notice.

Appendix

◆ Guidance and manufacturer's declaration – electromagnetic emissions

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	RF energy is used only to maintain device's operation. Therefore, its RF emissions are so low that it's not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not Applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable	

◆ Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	In the case of air discharge testing, the climatic conditions shall be within the following ranges: Ambient Temperature: 15°C–35°C, Relative Humidity: 30%–60%.
Power frequency magnetic field IEC 61000-4-8	30 A/m 50 or 60 Hz	30 A/m 50 or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

◆ Recommended separation distances between portable and mobile RF communication equipment and the device.

The device is intended for use in an electromagnetic environment where radiated RF disturbances are under control. User can help prevent electromagnetic interference by keeping the device at a minimum distance from portable and mobile RF communications equipment (transmitters). Below table details the maximum output power of transmitter:

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz Not Applicable	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3 \sqrt{P}$
0.01	Not Applicable	0.12	0.23
0.1	Not Applicable	0.38	0.73
1	Not Applicable	1.2	2.3
10	Not Applicable	3.8	7.3
100	Not Applicable	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

◆ **Guidance and manufacturer's declaration – electromagnetic immunity**

The device is intended for use in the electromagnetic environments listed below, and should only be used in such environments:

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3V rms At 0.15-80 MHz 6V rms At ISM & Radio Amateur Freq	Not Applicable	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	10 V/m at 80-2700 MHz	10 V/m at 80-2700 MHz	<p>Recommended separation distance</p> $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2.7 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>
Proximity fields from RF wireless communications equipment IEC 61000-4-3	AM Modulation And 9-28V/m at 385-6000MHz, Pulse Mode and other Modulation. The system shall be tested as specified in IEC60601-1-2 Table 9 for proximity fields from RF wireless communications equipment using the test methods specified in IEC 61000-4-3	AM Modulation And 9-28V/m at 385-6000MHz, Pulse Mode and other Modulation. The system shall be tested as specified in IEC60601-1-2 Table 9 for proximity fields from RF wireless communications equipment using the test methods specified in IEC 61000-4-3	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

To prove the accuracy of CTD504, we perform clinical experiments at one hospital by professionals on 3 groups of test subjects: new born to 1-year-old for Group A, 1-year-old to 5-year-old for Group B, and 5-year-old or above for Group C.

Each group consists of 50 effective test subjects, and at least 30 % of them are febrile, with the exception for those who are under 3-month-old.

All clinical test procedures are performed in compliance with EN 12470-5 and ISO 80601-2-56.

Statistic of subjects for Infrared Ear Thermometer

EN12470-5			
Statistic of Subjects	Repeatability(≤ 0.3)	Mean bias($\leq \pm 0.2$)	SD(≤ 0.4)
Group A	0.109	-0.088	0.232
Group B	0.104	0.017	0.267
Group C	0.104	-0.036	0.215

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- CITIZEN is a registered trademark of Citizen Watch Co., Ltd. Japan.
CITIZEN SYSTEMS JAPAN CO., LTD.