INSTRUCTION MANUAL FOR
COMPRESSOR NEBULIZER

GC806

FIG. A





Figure Relationship between the Cumulative Underize (%) and Particle Size

Donawa Lifescience
Piazza Albania, 10, 00153 Rome, Italy
<https://www.donawa.com/wli/main/contatti.index>

GLOBALCARE MEDICAL TECHNOLOGY CO., LTD
7th Building, 39 Middle Industrial Main Road, European Industrial
Zone, Xiaolan Town, 528415 Zhongshan City, Guangdong Province,
PEOPLE'S REPUBLIC OF CHINA
Phone : +86 760 22589901

CITIZEN is a registered trademark of Citizen Watch Co., Ltd. Japan.
Design and specifications are subject to change without notice.

English	GC806 COMPONENTS - (PIC. A)	
1. Compressor housing	6. Nebulizer holder	11. Adult mask
2. ON/OFF switch	7. Air-vent openings	12. Pediatric mask
3. Power cord	8. Nebulizer	13. Spare air filters
4. Air outlet	9. Mouthpiece	14. Nosepiece
5. Filter compartment	10. Air tube	15. Applied parts case

IMPORTANT: READ CAREFULLY THE INSTRUCTIONS FOR USE BEFORE OPERATION.

GENERAL INFORMATION

Intended purpose:

The device is an aerosoltherapy system suitable for domestic use.
The device is designed for the production of compressed air to operate a small volume nebulizer for the production of medical aerosol for lung respiratory disorders.

Patient population:

The device is intended for use with children from 2 years old, adolescent and adult patients.

Intended users:

The user is an adult person without specific knowledge or professional ability.
The patient is the intended operator except in case of children and patients that requires special assistance.

Indications:

Lung acute or chronic diseases of respiratory disorders organs, or inflammation of the upper respiratory system.

Contraindications:

The device is not indicated to use with quick-relief medications during life-threatening asthma attack. No contraindications exist to the administration of aerosols by inhalation. Contraindications related to the medicine used must be checked on the medicine package leaflet. Consult your physician in case of doubts.

Warnings:

- Use the device only as described in this manual and therefore as an aerosoltherapy system, following the indications of your doctor. Any use different from the intended one is to be considered improper and hence dangerous; the manufacturer cannot be held liable for any damage caused by improper, incorrect and/or unreasonable use, or if the equipment is connected to electrical installations which do not comply with current safety regulations.
- Keep this manual for future reference.
- Do not operate the unit in presence of any anaesthetic mixture inflammable with oxygen or nitrogen protoxide.
- The correct functioning of the equipment can be affected by electromagnetic interferences which exceed the limits indicated by the European standards in force. In case this device interferes with other electrical devices, move it and plug it to a different power socket.
- In case of failure and/or malfunction, read the "POSSIBLE PROBLEMS AND HOW TO SOLVE THEM" section. Do not tamper or open the compressor housing.
For repair operations address only to a technical service centre authorized by the manufacturer and require the use of original spare parts. The non-observation of the above mentioned indications can compromise the device safety.
- Use the dedicated spare parts only.
- Comply with the safety regulations concerning the electrical devices and in particular:
 - use only original components;
 - never submerge the unit in water;
 - never wet the device it is not protected against sprinkles;
 - never touch the unit with wet or moist hands;
 - do not leave the unit exposed to the weather elements;
 - place the unit on a stable and horizontal surface during its operation;
 - the use of this device by children and disabled requires always the close supervision by an adult with full mental faculties;
 - do not pull the power cord or the device itself to unplug it from the power socket;
 - the power plug in the separation element from the grid power; keep the plug accessible when the device is in use;
- Before plugging in the device, make sure that the electrical rating, shown on the rating plate on the bottom of the unit, corresponds to the mains rating.
- In case the power plug provided with the device does not fit your wall socket, address to qualified personnel for the replacement of the plug with a suitable one. In general, the use of adapters, simple or multiple, and/or extension cable is not recommended. If their use is indispensable, it is necessary to use that complying with safety regulations, paying however attention they do not exceed the maximum power limits, indicated on adapters and extension cables.
- Do not leave the unit plugged in when not in use; unplug the device from the wall socket when it is not operated.
- The installation must be carried out according to the instructions of the manufacturer. An improper installation can cause damage to persons, animals or things, for which the manufacturer cannot be held responsible.
- The power cord of this device cannot be replaced by the user. In case of a power cord damage, address to a technical service centre authorized by the manufacturer for its replacement.
- The power supply cord should always be fully unwound in order to prevent dangerous overheating.
- Before performing any maintenance or cleaning operation, turn off the device and disconnect the plug from the main supply.
- Some parts of the unit are so small that they may be swallowed by children; keep the equipment out from children's reach.
- If you decide not to use the device any longer, it is recommended to dispose of it according to the current regulations.
- Make sure to:
 - use this device only with medicines prescribed by your doctor.
 - make the treatment using only the applied part recommended by your doctor depending on the pathology.
 - use the nasal fork only if expressly indicated by your doctor and paying attention to NEVER introduce the bifurcations in the nose, but only bringing them as close as possible.
- Check in the medicine package leaflet for possible contraindications for use with common aerosoltherapy systems.
- Applied parts are single patient reuse. Device is multi-patient reuse.
- Nebulizing system is not suitable for use in anaesthetic breathing system or a ventilator breathing system. The patient cannot undergo MRI scan while using this device.
- Device is designed to nebulize solution and suspension liquids.
- Do not use the device, if you spot damage or you notice something unusual.
- This device consists of delicate components and must be handled with care. Observe the storage and operating conditions in the chapter "TECHNICAL SPECIFICATIONS".
- The use of this device is not a replacement for visiting the doctor. In case of no health improvement after treatment consult again your doctor.
- To avoid strangulation and entanglement, keep cable and air tubes out of reach of young children.

INSTRUCTIONS FOR USE

- The device must be checked before each use, in order to detect possible functioning anomalies and/or damages due to transport and/or storage.
- During inhalation, sit upright and relaxed at a table and not in an armchair, in order to avoid compressing your respiratory airways and impairing the treatment effectiveness.
- The applied parts must be used only with a single patient; it is not recommended to use them with several patients.
- After unpacking the device, check it for visible damages or defects; pay particular attention to cracks in the plastic housing, which may expose electrical components. Check for parts integrity.
- Before using the device, proceed with the cleaning operations as described in the "CLEANING AND DISINFECTION" section.
- Open the nebulizer by turning counterclockwise the top (Pic. B).
- Make sure that the medicine conduction cone is properly fitted on the air conduction cone inside the nebulizer (Pic. C).
- Put the prescribed quantity of medicine into the nebulizer (Pic. D).
- Close the nebulizer by turning clockwise the two parts, paying attention that they are thoroughly sealed (Pic. E).
- Connect one end of the air tube to the nebulizer (Pic. F) and the other end to the air outlet (Pic. G) on the device.
- If you need to use the mask, insert it directly onto the nebulizer (Pic. H).
- Plug the device into the wall socket, making sure that the main supply corresponds to the electrical rating of the device.
- To start the treatment, set the I/O switch into the « I » position.
- ATTENTION:** This device is for intermittent use: 30 min ON / 30 min OFF. Do not use the product continuously for 30 minutes or longer.
- Inhale the aerosol solution using the prescribed applied part.
- When the treatment has been completed, switch the unit off by setting the I/O into the « O » position and disconnect the plug from the wall socket.
- Wash the applied part as described in the "CLEANING AND DISINFECTION" section.
- The device requires no calibration.
- Tampering with the device is absolutely forbidden.
- No modification to the device is permitted.

CLEANING AND DISINFECTION

Device(s): Aerosoltherapy applied parts	
WARNINGS	Before the first use and after each treatment Follow carefully the cleaning and disinfection instructions of the applied parts as they are very important for the device performances and the therapy success. Use original parts only. Do not clean or disinfect the air tube. The applied parts can not be cleaned and disinfected by automated method. Do not boil nor autoclave the masks.
Limitations on processing	The nebulizer must be replaced after a long period of inactivity, in case it shows deformations or breakings, or when the nebulizer nozzle is obstructed by dry medicine, dust, etc. We recommend replacing the nebulizer after a period between 6 months and 1 year, depending on the usage. The maximum times of cleaning and disinfection of the nebulizer are 360 times. The maximum times of cleaning and disinfection of the mask, the nosepiece and the mouthpiece are 360 times.

INSTRUCTIONS

Preparation before cleaning	Detach the air tube from the nebulizer. Detach the mouthpiece or nosepiece or mask from the nebulizer. Make sure that all residual volume is removed from the nebulizer. Disassemble the nebulizer by turning counterclockwise the top and remove the medicine conduction cone. Note: To avoid microbial growth and residual drug drying, clean and disinfect the applied parts immediately after each use.
Cleaning:	Briefly rinse all parts in running tap water for at least 10 seconds beforehand. Mix little dish detergent and warm tap water (i.e. FAIRY brand dish washing liquid for hand washing in a ratio of 2ml: 1L) in a clean container. Immerse the components of disassembled nebulizer, the mask, the mouthpiece and the nosepiece into the mixed water for about 5 min. And then scrub all the surface of all the components with clean and small brush for at least 8 times. Afterwards rinse all parts thoroughly in running tap water at least 30 seconds, to fully remove any possible dish washing liquid residue.
Disinfection	After cleaning, disinfect all of the disassembled parts (only parts that have been cleaned can be disinfected effectively). Boil the disassembled nebulizer, the mouthpiece and the nosepiece for 5 minutes in boiling tap water. Immerse the mask in a 2% sodium hypochlorite (NaOCl) solution for 15 minutes, (i.e. a solution made from the disinfecting agent Amuchina® or a 2% sodium hypochlorite solution mixed by your pharmacist). Afterwards immerse the mask with sterile water for 3 min, and then rinse all the surface of the mask with other clean sterile water twice, to fully remove any possible residue of the disinfectant solution.
Drying	Reassemble the nebulizer components and connect it to the air-outlet, switch the device on and let it work for 10-15 minutes. Let all parts dry completely before reassembling and reuse to avoid risk of germ grow. Care needs to be taken not to contaminate the parts after they have been cleaned and disinfected. Contamination can be avoided by good hand washing and not touching the inside sections of the device when laying them out to dry or when reassembling.
Inspection	Inspect all product components after each cleaning and disinfection. Replace any broken, misshapen or seriously discolored parts.
Packaging	Pack dry parts in a clean and sealed container when not in use. Do NOT pack wet or damp parts.
Storage	Storage conditions refer to "TECHNICAL SPECIFICATIONS"
Transportation	Note: re-clean and re-disinfect the parts if they are stored more than one day. After cleaning and disinfection, always transport the parts with clean and sealed container. Contamination can be avoided by good hand washing and not touching the inside sections of the parts when take out and re-assemble the parts for use.

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

△ FILTER REPLACEMENT

In normal usage conditions, the air filter must be replaced approximately after 200 working hours or once a year. We recommend to periodically check the air filter (10-12 treatments) and, if the filter shows a grey or brown colour or is wet, replace it. Remove the filter and replace it with a new one.

Do not attempt to clean the filter to reuse it. The air filter shall not be serviced or maintained while in use with a patient.

△ Only use original filters. Do not use the device without the filter.

REPLACEMENT OF THE APPLIED PARTS

No. 8 to 14 in the figure are consumables. They are available to be purchased separately (model number: AT-GC202). When purchasing any of the above items, contact your local distributor.

△ Use original filters only.

△ Do not use the device without filter.

POSSIBLE PROBLEMS AND HOW TO SOLVE THEM

PROBLEMS:

- The device does not switch on
- The device does not nebulize or nebulizes weakly

SOLUTIONS:

- Make sure the power plug is firmly fitted to the wall socket;
 - Make sure that the device has been operating within limits indicated in this manual (30 min ON / 30 min OFF).
 - Make sure that the ends of the air tube are tightly fitted to the main unit and nebulizer;
 - Verify that the nebulizer is not empty or has been filled with the proper quantity of medicine (max 8 ml).
 - Verify that the nebulizer nozzle is not obstructed.
- Should the device not start working properly again, address to local distributor.

MAINTENANCE AND REPAIRS

In case of failure, address to local distributor. Do not open the device in any case. The unit has no user-serviceable parts within and does not need internal maintenance or lubrication.

TECHNICAL SPECIFICATIONS

Item	Specification
Max Pressure (acceptable range)	200-320 kPa
Max Pressure (average value)	220 kPa
Operating Pressure (with calibrated nozzle, Ø 0,50)	80 ~ 155 kPa
Protection type IP	IP 21
Sound noise level (Distance: 1m, Front direction) max level	≤58 dBA
Sound noise level -average level	≤52 dBA
Equipment conditions of use	30 ON / 30 OFF min
Dimensions	130x130x117 mm
Weight	1.3±0.2 kg
Service life	1000 Hours
Power supply	GC806-EU AC230V 50Hz - GC806-ME AC230V 50Hz - GC806-SA AC220V 60Hz - GC806-MX AC120V 60Hz -
Current ABSORPTION	≤ 1000 mA
Neb Rate (average value)	0.35 ml/min
Neb Rate (acceptable range)	0.23~0.60 ml/min
Nebulizer kit	Clenny
AEROSOL PERFORMANCES ACCORDING TO EN ISO 27427:2019 based on adult ventilatory pattern with Sodium Fluoride (NaF):	
Aerosol output	0.11 ml
Aerosol output rate	0.07 ml/min
Percentage of fill volume emitted per min:	3.5 %
Residual Volume	1.68 ml
Particle size (MMAD)	3.07 µm
GSD (geometric standard deviation):	2.03
RF (respirable fraction <5µm):	75 %
Large particle range (>5µm):	25 %
Middle particle range (2 to 5µm):	47.8 %
Small particle range (<2µm):	27.2 %
Min & Max filling volume:	2 ~ 8 ml
Operating Temperature	10 ~ 40 °C
Operating Humidity	10 ~ 95 % RH
Operating Pressure	700 ~ 1060 hPa
Storage Temperature	25 ~ 70 °C
Storage Humidity	10 ~ 95 % RH
Storage Pressure	700 ~ 1060 hPa
Protection against electrical shock	2 /

• Class II device referred to protection against electric shock.

• Nebulizer, mouthpiece, masks, and nosepiece are type BF applied parts.

Note: Technical specifications subject to change without notice.

AEROSOL CHARACTERISTICS IN ACCORDANCE WITH REGULATION EN ISO 27427

Standards applied:

Electric safety standards EN 60601-1 Electromagnetic compatibility according to EN 60601-1-2

The device is a Class IIa medical device. **Device Complies with European Regulation on Medical Devices EU MDR 2017/745**

DEVICE SYMBOLS

REF	Reference number	UDI	Unique device identifier	☒	Humidity limitation
I/O	ON / OFF	EC REP	Authorized representative	☒	Temperature limitation
~	Alternating current	Manufacturer	Manufacturer	☒	Ambient pressure limitation
Class II	Class II	SN	Serial number (LOT/serial number)	☒	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive - WEEE
Type BF	Type BF	Date of manufacture (year, month and day)	Date of manufacture (year, month and day)	#	Type number
CE 0123	Device Complies with European Regulation on Medical Device EU MDR 2017/745	CAUTION	CAUTION	MD	Medical device
IP21	Protection against solid foreign objects and against harmful effects due to the ingress of water.	Follow instruction for use	Follow instruction for use	☒	Importer
Single patient multiple use (for applied parts only)	Single patient multiple use (for applied parts only)	Corrugated cardboard recycling symbol	Corrugated cardboard recycling symbol		

△ WEEE- Waste Electrical and Electronic Equipment

This product is subject to European Directive 2012/19/EU on waste electrical and electronic equipment and is marked accordingly. Never dispose of electronic devices with household waste. Please seek out information about the local regulations with regard to the correct disposal of electrical and electronic products. Correct disposal helps to protect the environment and human health. Please report any serious incident that has occurred in relation to the device injury or adverse event to the local competent authority and to the Manufacturer or to the European Authorised Representative (EC REP)
Vigilance contacts points: <https://ec.europa.eu/growth/sectors/medical-devices/contacts/>

DISPOSAL

Item	Component	Disposal Description	Photo
1	Device	The component mainly includes plastic and electronic components. All comply with RoHS and REACH, and all could be safely disposal. This product is subject to European Directive 2012/19/EU on waste electrical and electronic equipment and is marked accordingly. Never dispose of electronic devices with household waste. Please seek out information about the local regulations with regard to the correct disposal of electrical and electronic products.	
2	Applied parts (air tube, nebulizer, nosepiece, mouthpiece, masks)	The components are plastic. All comply with RoHS and REACH. Before disposal, all the parts should be cleaned according to manual instructions and then disinfected by boil for 5 minutes.	

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