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

Device Labeling

For Digital Thermometer

Model: CT422 (With Buzzer & Water Resistance)

Please read this manual before use

Parts Names

Probe	(Temperature measuring part)
Display	(Display Screen)
Push Button	(On/Off Switch)
Probe tip (Fever sensor)	(Heat/Fever sensor)
Main Body	(Using antibacterial material)
Battery Cover	(Part that covers battery)

Precautions

1. Do not use this thermometer for measurements other than body temperature. The thermometer will be damaged if it is used to measure the temperature of baths, refrigerators, etc.
2. Do not drop the thermometer and do not expose the thermometer to heavy shock because the thermometer is not designed with a shock resistant system.
3. Do not bend the thermometer.
4. Do not bite the thermometer.
5. Do not disassemble the unit, do not make any modifications and do not touch the circuitry.
6. Do not spill chemical solvents on the unit.
7. Do not store the unit in direct sunlight, under high temperature, or high degree of humidity or dust.
8. Do not boil the probe.
9. Do not use the thermometer if it operates irregularly or shows signs of malfunction.
10. Do not use thinner or benzene to clean the thermometer.
11. When using the thermometer after storage at low temperature, place it in a warm place for 5 to 10 minutes before use.
12. Keep the thermometer out of reach of children.
13. Wipe the thermometer clean before storing.
14. For best use, it is recommended to allow over 5 minutes measuring time for oral temperatures (mouth) and over 10 minutes for underarm (axillary) temperatures.

Important Data

Instruction	: Use to measure body temperature
Measurement Range	: 32°C – 42.9°C
Accuracy	: ±0.1°C
Display	: 3 digits with decimal point, M (Memory mark), °C (°C indicator), 🔋 (Battery symbol)
Rated voltage	: DC 1.5V (Direct current)
Battery life	: Approximately 3 years (if used 1 time a day, for 10 minutes each time)
Weight	: Approximately 10g (including battery)
Operating environment	: Temperature 10 to 40°C, Relative humidity 30 to 85%
Storage/transport environment	: Temperature -20°C to 60°C, Relative humidity 95% or less Internal power unit
Electrical shock protection	: Internal power unit

Medical Device Documentation Date: 20/03/2018

How to Operate

1. Press the Power Button to turn on the device. The screen will display **88.8°C** to test its readiness for effective use. If you press and hold the Power Button, the screen will display the last measured value, **36.8°C**.
2. When you release the Power Button, **36.5°C** will appear on the screen, and then **L o °C** and °C will flash, indicating that the thermometer is ready for measurement.
3. When the temperature stabilizes, °C will stop flashing and an alarm will sound.
4. Press the Power Button to turn off the device, or the device will automatically turn off after 30 minutes.

Notes

1. If the ambient temperature is higher than 32°C, the display will not show **L o °C** after showing **36.5°C**.
2. If the temperature is lower than 32°C, the display will only show the symbol "**L o °C**". However, if the temperature is higher than 42.9°C, only the symbol "**H, °C**" will be displayed.
3. Temperature measurements should be taken under adult supervision, and rectal temperature measurements with this thermometer are not recommended.

Temperature Measurement

1. **Oral Temperature Measurement**
Place the thermal sensor under the tongue and close your mouth tightly for an accurate reading. Generally, oral temperatures are about 0.2–0.5 degrees Celsius higher than under the armpit.
2. **Underarm Temperature Measurement**
Place the device snugly under the armpit, ensuring the temperature sensor is in direct contact with the skin. Place your hand on your chest to keep your arm securely against the device. Note: The resulting temperature is usually about 0.5 degrees Celsius lower than the oral reading. Caution: Always use the thermometer to measure in the same location and at the same time each day.

Battery

When the battery is low, a low battery indicator  will appear and flash on the display. Replace the battery.

Replacing the Battery

1. Use a screwdriver to slide open the battery compartment cover by clicking the hook located  on the back of the device.
2. Pull out the old battery and insert the new one (+).
3. Ensure the positive (+) side is visible for correct positioning. Close the battery compartment cover in the same direction.

Caution:

- Swallowing batteries is life-threatening. If you encounter any problems, seek immediate medical attention.
- Do not throw batteries into a fire as they may explode.
- Do not recharge, disassemble, or short-circuit old batteries.

Cleaning the instrument

Clean the thermometer using a cloth dampened with ethyl alcohol or a neutral cleaning solution. Do not use thinner or gasoline.

For the EMD information, please scan the 2D code

Importer : Eminence International Co., Ltd.,
5 Prachanukool 3 Alley (Ratchadaphisek 66),
Ratchadaphisek Road, Wongsawang Subdistrict,
Bang Sue District, Bangkok 10800, Thailand.
Tel: 02 910 1255.
Medical device importer, registration number: SN. 742/2554.

Manufacturer : CITIZEN SYSTEMS (JIANGMEN) CO., LTD.,
Building 6, 399 Jinxing Road, Jianghai District, Jiangmen,
Guangdong, China.

Product Owner : CITIZEN SYSTEMS JAPAN CO. LTD., Japan.