INSTRUCTION MANUAL FOR DIGITAL BLOOD PRESSURE MONITOR

CH-650

English
Français
中国語
Русский
Español
Português
한국어

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Contents

- 3 - GENERAL REMARKS
- 5 - PRECAUTIONS FOR USE AND MAINTENANCE
- 10 - IDENTIFICATION OF PARTS
- 12 - LOADING THE BATTERIES
- 13 - REMOVING THE BATTERIES
- 14 - HOW TO WRAP THE WRIST CUFF
- 16 - TIPS FOR OBTAINING ACCURATE MEASUREMENTS
- 18 - MEASURING YOUR BLOOD PRESSURE
- 20 - CALLING UP THE DATA STORED IN MEMORY
- 22 - ERROR INDICATIONS
- 23 - ABOUT "BLOOD PRESSURE"
- 24 - BLOOD PRESSURE Q&A
- 25 - BEFORE REQUESTING REPAIRS OR TESTING
- 26 - SPECIFICATIONS

Ensure that you have all the following components


Blood pressure
monitor unit

Storage case

(for monitor)
2 AAA batteries

Instruction
Manual

GENERAL REMARKS

General remarks on blood pressure and blood pressure measurement


1. To deflate the cuff quickly, press the "START/STOP" switch.
2. The attached cuff is suitable for wrist circumferences of 13.5-21.5 cm.
3. When attaching the cuff, the body of a blood pressure monitor should be placed on the inside of your wrist.
4. The monitor is not intended to be a diagnostic device, because it is intended for home use.
5. Self measurement is not therapy! In no event should you change the dosage of your medication prescribed by your physician.
6. As preparation for measuring your blood pressure, you should urinate and then remain relatively still for 10 to 15 minutes prior to measurement.
7. Exercising, eating, drinking, smoking, etc., prior to measurement may affect the results.
8. Your blood pressure varies constantly throughout the day. Measurements should be taken regularly at the same time each day.
9. DO NOT place too much emphasis on the results of one measurement. Keep a continuing record of your blood pressure variations. A complete picture can only be obtained from a large number of readings.
10. Emotional stress tends to cause your blood pressure to rise.

To obtain correct blood pressure measurements

1. Take five or six deep breaths and then relax before measuring your blood pressure. If you are tense when taking the measurement, you will not get a valid reading.
2. Your blood pressure will be elevated if you are anxious or irritated, suffering from lack of sleep or constipation, or have just taken some exercise or eaten a meal.
3. Wrap the cuff around your wrist snugly. ( See page 14)
4. DO NOT measure your blood pressure after bathing or drinking.
5. If you feel the need to urinate, do so before measuring your blood pressure.
6. Measure your blood pressure where the room temperature is around 20°C. DO NOT measure your blood pressure when it is too cold (below 10°C) or too hot (above 40°C) in the room.
7. DO NOT try to measure your blood pressure immediately after drinking coffee or tea or after smoking.
8. Measure your blood pressure when you are relaxed and still. Keep the center of the cuff at the level of your heart and DO NOT move your arm or talk.

GENERAL REMARKS

9. DO NOT measure blood pressure when the cuff has been on your wrist for a few minutes or more. During this time, your wrist will have built up a higher blood pressure and you will not get valid reading.
10. If the Body Movement Indicator or Irregular Heartbeat (IHB) Indicator is displayed, relax and measure your blood pressure again. If the Irregular Heartbeat (IHB) Indicator appears frequently, you are advised to consult a doctor about your health condition. (👉 See page 19)
11. The figure below shows the blood pressure classifications under the WHO standards.

Remarks:

The graph is not exact, but may be used as a guide in understanding non-invasive blood pressure measurements. The device is only intended for use by adults.

Definitions and Classifications of Blood Pressure Levels

12. Intended Use
 - This device is noninvasive blood pressure monitor by oscillometric method. This device can measure the systolic blood pressure (SYS), the diastolic blood pressure (DIA) and the pulse rate automatically.
 - This device is intended to be used for checking personal health condition at home under the direction of a physician and is not intended to be a diagnostic device.
 - This device is not intended for use with neonates or infant, and not intended for automatic cycling measurement.

Measure your blood pressure at the same time each day.

- * Your blood pressure changes all the time. This means that data gathered over a long period has far more significance than data from just one measurement. For this reason, you must measure your blood pressure on a daily basis. Ideally, you should measure your blood pressure at the same time each day, wherever possible.

PRECAUTIONS FOR USE AND MAINTENANCE

Precautions for use

1. If you suffer from heart disease, high blood pressure or other circulatory disease, consult your physician before using the monitor.
2. If you feel abnormal pressure of the cuff during use or any other irregularity, reduce the pressure by immediately switching off the power or remove the batteries and consult the sales outlet where purchased.
3. If you think the measurement is abnormal or if measurement makes you feel unwell, discontinue use and consult your physician.
4. Blood pressure measurement may not be possible for anyone with a weak pulse or arrhythmia.
5. Repeated blood pressure measurement may cause problems such as congestion or swelling in some people.
6. Frequently repeated blood pressure measurements will not give accurate results. Allow an interval of about 1 minute between measurements.
7. If you suffer from a severe problem with blood circulation in your arms, consult your physician before using the monitor. Failure to do so could be hazardous to your health.
8. Measurement may not be possible for anyone with insufficient blood flow to the area where measurements will be taken or who suffers from a frequent irregular heartbeat. Consult your physician for advice on whether to use the monitor.
9. DO NOT wrap the cuff around an injured wrist.
10. DO NOT wrap the cuff around a wrist in which a drip (intravenous infusion) is inserted, or which is being used for blood transfusion as part of medical treatment. Doing so could result in an injury or a serious accident.
11. DO NOT use the monitor in the vicinity of flammable gases such as those used for anaesthesia. Doing so could ignite the gases and cause an explosion.
12. DO NOT use the monitor in enriched oxygen environments such as a hospital's hyperbaric chamber or oxygen tent. Doing so could ignite the oxygen and cause a fire.

PRECAUTIONS FOR USE AND MAINTENANCE

13. DO NOT use mobile phones near the monitor as this could result in a malfunction.
14. If you use a cardiac pacemaker, consult your physician before using the monitor.
15. Be sure to use this unit only for measuring blood pressure. DO NOT use it for any other purpose.
16. Be sure to use only our authorized parts and accessories. Parts and accessories not approved for use with the device may damage the unit.
17. DO NOT use this unit on infants.
18. Blood pressure measurement may not be possible for anyone with common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation.

Maintenance Precautions

1. DO NOT store the blood pressure monitor in locations exposed to direct sunlight, high temperatures (over 60°C), low temperatures (below -20°C), high relative humidity (over 95%) or excessive amounts of dust.
Make sure to store the blood pressure monitor, where children, pets and or pests are not there.
2. DO NOT drop the blood pressure monitor or subject it to other shocks or vibration.
3. Remove the batteries if the monitor will be left unused for a long period.
4. DO NOT attempt to disassemble the monitor.
5. DO NOT bend the cuff excessively.
6. If the monitor is very dirty, wipe it clean with a cloth moistened with sterilizing alcohol or a neutral detergent. Then wipe it with a dry cloth.
7. NEVER clean the blood pressure monitor with thinners or benzene, as this could damage the monitor.
8. To clean the cuff, wipe it with a moist cloth. Avoid hard rubbing, as this will cause air leakages.

PRECAUTIONS FOR USE AND MAINTENANCE

Warning

No modification of this device is allowed.

Symbols Explanation

: The CE marking is meant to provide information to market inspectors in the EU member countries.

: Type BF applied part

: Keep dry

: Warning

: Refer to instruction manual before use

: Appliance compliance WEEE directive

: Manufacturer

Warning

The device has been tested and homologated in accordance with EN60601-1-2 for EMC. This does not guarantee in any way that the device will not be affected by electromagnetic interference. Avoid using the device in high electromagnetic environment.

PRECAUTIONS FOR USE AND MAINTENANCE

Nature and frequency of maintenance:

This product is designed for use over an extended period of time; however, it is generally recommended that it be inspected every two years to ensure proper function and performance.

Protect the nature environment:

Please help to protect natural environment by respecting national and/or local recycling regulations when disposing of the battery and the product at the end of their useful live.

WEEE MARK

If you want to dispose this product, do not mix with general household waste. There is a separate collection systems for used electronics products in accordance with legislation under the WEEE Directive (Directive 2012/19/EU) and is effective only within European Union.

Information on Disposal in other Countries outside the European Union.

If you wish to discard used batteries, please contact your local authorities or dealer and ask for the correct method of disposal.

Note for the battery symbol.

The symbol might be used in combination with a chemical symbol. In this case it complies with the requirement set by the Directive for the chemical involved.

Reference European standard:

The blood pressure measuring device corresponds to regulation EN60601-1, EN1060-1, EN1060-3, EN1060-4.

Calibration:

The blood pressure measuring device is generally recommended to have the monitor inspected every two years to ensure correct functioning and accuracy. Please contact a distributor.

PRECAUTIONS FOR USE AND MAINTENANCE

Name : CITIZEN SYSTEMS JAPAN CO.,LTD.
Address : 6-1-12, Tanashi-cho, Nishi-Tokyo-shi, Tokyo 188-8511, Japan

Factory

Name : CITIZEN SYSTEMS (JIANGMEN) CO., LTD.
Address : Building 6, 399 Jinxing Road, Jianghai District, Jiangmen, Guangdong, China

European Representative:

Name : EMERGO EUROPE
Address : Prinsessegracht 20, 2514 AP The Hague, The Netherlands

- For technical data, please contact our European Representative

Made in China


IDENTIFICATION OF PARTS

Component names

<Main Unit>



<Display>

LOADING THE BATTERIES

(Batteries supplied with the unit should be loaded in compartment before you use the blood pressure monitor.)

1 Open the cover of battery compartment.


Pressing down on the battery cover, slide the cover off in the direction of the arrow.

2 Place batteries in the compartment.


Put the batteries paying attention to the positive and negative terminal symbols \oplus and \ominus .

* The battery goes in negative side \ominus first.


3 Close the battery cover.

Slide the cover in the direction of the arrow. Push the cover on until it clicks.

REMOVING THE BATTERIES

- When the icon  is displayed or nothing appears in the display, replace the old batteries with two new ones. All two batteries should be changed at once.
- Use only alkaline batteries.
- Do not use rechargeable batteries.
- Do not use alkaline (LR03) and manganese (R03) batteries together.
- Batteries included are for demonstration purposes only. Battery life may be shorter than specified.
- When disposing of used batteries, comply with governmental regulations or environmental public institution's rules that apply in your country/area.

HOW TO WRAP THE WRIST CUFF


Wrap the wrist cuff around your wrist properly. Improper wrapping of the cuff may result in measuring errors.

- 1 Open the wrist cuff.**
- 2 Wrap the cuff around your left wrist so that the monitor is positioned on the inside of your wrist.**

- Wrap the wrist cuff on bare skin.
- The wrist cuff is designed to fit wrists from 13.5 to 21.5 cm in circumference.


- 3 Loosely fit the wrist cuff by pulling the end of the cuff.**


HOW TO WRAP THE WRIST CUFF

- 4** Position the top edge of the cuff approximately 1 to 2 cm away from the crease of your wrist.

- 5** Secure the wrist cuff with the hook-and-loop fastener.

- If the wrist cuff is wrapped too loose or tight, errors or incorrect readings may result.

You may use your right wrist to measure your blood pressure.

However, your blood pressure may differ by about 10mmHg depending on whether measured on your left or your right wrist, so ideally you should measure your blood pressure on the same wrist each time.

To obtain accurate measurements

- Wrapping the cuff too loosely or too tightly may prevent correct measurement.
- Sit up straight while taking the measurement.

TIPS FOR OBTAINING ACCURATE MEASUREMENTS

Your blood pressure varies according to your posture, the time of day and a range of other factors. Ideally, you should measure your blood pressure in the same posture at the same time every day.

- ① Make sure that your posture is correct.
Sit up straight, take a deep breath and relax.

- ② Place your palm facing up and keep your fingers relaxed.

- ③ It is important to have the cuff level with your heart.

Place a folded towel or similar under your arm.

TIPS FOR OBTAINING ACCURATE MEASUREMENTS


To obtain correct measurements

- Do not clench your fist tightly.

- Do not wrap the wrist cuff too loosely.

- Do not change the position of the monitor.


- Do not move your body or talk during measurement

MEASURING YOUR BLOOD PRESSURE


- 1 Wrap the wrist cuff around your left wrist.**
( See page 14 "HOW TO WRAP THE WRIST CUFF".)
- 2 Make sure that your posture is correct.**
( See page 16 "TIPS FOR OBTAINING ACCURATE MEASUREMENTS".)
- 3 Press the "START/STOP" switch to start measurement.**

- All digits displayed is an initial display, indicating that the monitor is functioning normal.
The  mark will also be displayed but this does not indicate that the battery needs to be replaced.

 is displayed when a pulse is detected.

4 The measurement results are displayed.

The Body Movement Indicator is displayed on the measurement results display if you move your hand or wrist during measurement, causing a substantial pressure change to be detected.

- If the Body Movement Indicator is displayed, measure your blood pressure again.

If an irregular heartbeat is detected during measurement, the Irregular Heartbeat (IHB) Indicator is displayed after measurement is completed.

An irregular heart beat is defined as a heart beat that varies by $\pm 25\%$ from the middle of 5 heart beat intervals during the blood pressure measurement.

- Measurement may not be correct if your heartbeat fluctuates greatly during measurement. If the Irregular Heartbeat (IHB) Indicator is displayed, measure your blood pressure again.

Caution

If the Irregular Heartbeat (IHB) Indicator appears frequently, you are advised to consult a doctor about your health condition.

MEASURING YOUR BLOOD PRESSURE

5 Finishing measurement.

- Press the “START/STOP” switch to turn the monitor off.

AUTOMATIC POWER OFF FEATURE

- After taking the blood pressure measurement, pressing the “START/STOP” switch will turn the monitor off. However, if you do not press the switch, the monitor’s Auto OFF function will automatically turn the monitor off after 3 minutes.


CALLING UP THE DATA STORED IN MEMORY

60 measurement results can be stored in the memory.

The average value is calculated automatically to help you manage your daily health.

HOW TO CALL UP THE DATA

1 Press the “MEMORY” switch.

- “3A” indicating the average value appears on the left section of the display. Then, the average value, based on the last three measurements, is displayed. (“3A” and the average value do not appear if there are two or fewer measurements stored in memory.)


ERROR INDICATIONS

DISPLAY MARK	CONDITION/CAUSE	CORRECTIVE ACTION
	Appears when the wrist cuff is wrapped too loose.	Re-wrap the wrist cuff snugly. Take another reading, making sure you remain still. (👉 See page 14 "HOW TO WRAP THE WRIST CUFF".)
	Appears when your wrist moved during the measurement.; or Appears when the cuff is not wrapped correctly.	Re-wrap the wrist cuff correctly. Take another reading, making sure you remain perfectly still. (👉 See page 14 "HOW TO WRAP THE WRIST CUFF".)
	Appears when the batteries have run out.	Replace two alkaline batteries with new ones. (👉 See page 13 "REMOVING THE BATTERIES".)
	Appears when the pulse rate falls outside the measuring range (below 39 beats/min. or above 181 beats/min.)	Re-wrap the wrist cuff properly. Take deep breaths to relax, and measure again. (👉 See page 16 "TIPS FOR OBTAINING ACCURATE MEASUREMENT".)
	The unit does not function properly.	Contact the store where it was purchased or your local service center. (👉 See page 25 "BEFORE REQUESTING REPAIRS OR TESTING".)

ABOUT "BLOOD PRESSURE"

What is blood pressure?


The heart is a pump that circulates blood throughout your body. Blood is pumped from the heart at a constant pressure into arteries. This pressure is called the arterial blood pressure and represents, in general terms, your blood pressure. Blood pressure is indicated by several kinds of pressures, including the systolic pressure that occurs when the heart pumps blood, and the diastolic pressure that occurs when blood returns to the heart.

Your blood pressure changes all the time.

Your blood pressure differs according to your age, gender and a range of other factors. It is likely to be affected by your biorhythm during the day and by your posture, physical activities, mental activities, level of stress and even by the air temperature. The blood pressure of a healthy person generally varies within a day.

[Example of blood pressure variations within a day]

(Blood pressure values measured at five-minute intervals through the day)

● Bevan AT, Honour AJ, Stott FH. Clin Sci 1969;36:329-44.

BLOOD PRESSURE Q&A

Q Why is the blood pressure measured at home different from that measured by the physician or at a hospital?

A Your blood pressure may vary due to exercise, the ambient temperature or your mental state. When you have a physician or nurse measure your blood pressure, the pressure value is likely to be 10 to 20 mmHg higher than usual due to anxiety and/or stress. Knowing what affects your blood pressure helps you to better monitor your health.

Q Why does the blood pressure obtained vary with measurement?

A Our blood pressure is adjusted by automatic nerve function. Blood pressure differs with every heart beat. We are likely to think that our blood pressure is constant, but it varies if you take measurements in repetition. Blood pressure is susceptible to the time of day, such as morning and afternoon, season and atmospheric temperature. In addition, blood pressure is easily affected by mental stress or emotional ups and downs. It tends to increase when you feel tension or decrease when you are relaxed.

Q What is the benefit of measuring my blood pressure at home?

A Blood pressure measurements taken at home give reliable data as they can be taken when you are in a stable and relaxed condition. Physicians place considerable importance on blood pressure measurements taken at home. You should measure your blood pressure at the same time each day and record the variations in your blood pressure.

Always consult your physician for an interpretation of your blood pressure readings and to determine the proper treatment.

BEFORE REQUESTING REPAIRS OR TESTING

Make sure to check the following before sending your blood pressure monitor away to be tested or repaired.

Problem	Checks	Response
Display stays blank when you press the "START/STOP" switch.	Check whether the batteries are dead.	Replace all the batteries with new ones.
	Check whether the batteries are installed the right way round (⊕ and ⊖ oriented correctly).	Load the batteries in the correct direction.
No measurements can be taken.	Check whether the measurement icon  is displayed.	Attach the cuff properly.
	Check whether the cuff is attached correctly.	
	Did you stay calm during measurement?	Measure your blood pressure again, making sure to remain still.
	Note that blood pressure measurement may not be possible for someone with an extremely weak pulse or cardiac arrhythmia (irregular pulse).	
Your blood pressure readings are abnormally high or low compared with those taken in hospital.		 See page 24.
Other phenomena		Remove the batteries from the unit and change them with new ones.

SPECIFICATIONS

Model Number		CH-650
Measurement system		Oscillometric method
Display		Digital display type
Wrist circumference range		13.5 to 21.5 cm
Measurement range	Pressure	0 to 280 mmHg
	Pulse	40 to 180 pulse/min
Accuracy	Pressure	±3 mmHg
	Pulse	±5% of reading
Inflation		Automatic inflation by internal pump
Rated voltage		3V DC $\overline{\text{---}}$ ($\overline{\text{---}}$: direct current)
Exhaust		Electromagnetic quick exhaust valve
Power supply		2 AAA batteries (LR03)
Battery duration		Approx. 400 times (180 mmHg, once/day, 22°C)
Operating conditions	Temperature	10°C to 40°C
	Humidity	15% to 85% RH
Storage conditions	Temperature	-20°C to 60°C
	Humidity	10% to 95%RH
Memory Function		Up to 60 measurements
Electric shock Protection		Internal power unit
Degree of Protection		\boxtimes Type BF applied part
Mode of operation		Continuous operation
Protection against ingress of water		IPX0
Dimensions		65(W) x 74(H) x 29(D) mm
Weight		Approx. 103g w/o batteries
Service Life		5 years
Accessories		Set includes 2 AAA batteries (LR03) for the monitor, Instruction Manual, Storage case

* Applied part for this device is Cuff.

* A range in barometric pressure 700 hPa to 1060 hPa

* Overall system accuracy of this device is met with the requirement of item 7.9 of EN1060-3.

EN1060-3 : Non-invasive sphygmomanometers Part 3 (European Standard)

ELECTROMAGNETIC COMPATIBILITY INFORMATION

- WARNING**
- Portable RF communications equipment should be used no closer than 30 cm (12 inches) to any part of the [CH-650], including cables specified.
 - Use of this equipment adjacent to or stacked with other equipment should be avoided.
 - Use of accessories and options other than those specified (other than CITIZEN original parts) could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment.

Guidance and manufacturer's declaration - electromagnetic emissions		
The [CH-650] is intended for use in the electromagnetic environment specified below. The customer or the user of the [CH-650] should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The [CH-650] uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The [CH-650] is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	N/A	

Guidance and manufacturer's declaration - electromagnetic immunity			
The [CH-650] is intended for use in the electromagnetic environment specified below. The customer or the user of the [CH-650] should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact	±8 kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
	±15 kV air	±15 kV air	
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	N/A	N/A
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to earth	N/A	N/A
Voltage dips, short interruptions and voltage variations on power supply IEC 61000-4-11	0% U_T 0.5 cycle	N/A	N/A
	0% U_T 1 cycle		
	70% U_T 25/30 cycle		
	0% U_T 250/300 cycle		
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: U_T is the A.C. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity

The [CH-650] is intended for use in the electromagnetic environment specified below. The customer or the user of the [CH-650] should assure that it is used in such an environment.

Conducted RF IEC 61000-4-6	3 Vrms, 6 Vrms 150 kHz to 80 MHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the [CH-650], including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance N/A
	Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^{a)} , should be less than the compliance level in each frequency range. ^{b)} Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a). Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the [CH-650] is used exceeds the applicable RF compliance level above, the [CH-650] should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the [CH-650].
- b). Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the [CH-650]

The [CH-650] is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the [CH-650] can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and [CH-650] as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)			
	150 kHz to 80 MHz N/A	150 kHz to 80 MHz N/A	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3 \sqrt{P}$
0.01	N/A	N/A	0.12	0.23
0.1	N/A	N/A	0.38	0.73
1	N/A	N/A	1.2	2.3
10	N/A	N/A	3.8	7.3
100	N/A	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration - electromagnetic immunity

The [CH-650] is intended for use in the electromagnetic environment specified below. The customer or the user of the [CH-650] should assure that it is used in such an environment.

Immunity to proximity fields from RF wireless communications equipment IEC 61000-4-3	Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
	385	380-390	TETRA 400	Pulse modulation b) 18 Hz	1.8	0.3	27
	450	430-470	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710	704-787	LTE Band 13, 17	Pulse modulation b) 217 Hz	0.2	0.3	9
	745						
	780						
	810	800-960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	2	0.3	28
	870						
	930						
	930						
	1720	1700-1900	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0.3	28
	5240	5100-5800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0.2	0.3	9
	5500						
	5785						

NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the [CH-650] may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

- a). For some services, only the uplink frequencies are included.
 b). The carrier shall be modulated using a 50 % duty cycle square wave signal.
 c). As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

- CITIZEN is a registered trademark of Citizen Watch Co., Ltd. Japan.
 - Design and specifications are subject to change without notice.
-

CITIZEN SYSTEMS JAPAN CO., LTD.

6-1-12, Tanashi-cho, Nishi-Tokyo-shi,
Tokyo 188-8511, Japan

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SPECIFICATIONS

Model Number		CH-650
Measurement system		Oscillometric method
Display		Digital display type
Wrist circumference range		13.5 to 21.5 cm
Measurement range	Pressure	0 to 280 mmHg
	Pulse	40 to 180 pulse/min
Accuracy	Pressure	± 3 mmHg
	Pulse	± 5% of reading
Inflation		Automatic inflation by internal pump
Rated voltage		3V DC  ( : direct current)
Exhaust		Electromagnetic quick exhaust valve
Power supply		2 AAA batteries (LR03)
Battery duration		Approx. 400 times (180 mmHg, once/day, 22°C)
Operating conditions	Temperature	10°C to 40°C
	Humidity	15% to 85% RH
Storage conditions	Temperature	-20°C to 60°C
	Humidity	10% to 95%RH
Memory Function		Up to 60 measurements
Electric shock Protection		Internal power unit
Degree of Protection		 Type BF applied part
Mode of operation		Continuous operation
Protection against ingress of water		IPX0
Dimensions		65(W) x 74(H) x 29(D) mm
Weight		Approx. 103 g w/o batteries
Service Life		5 years
Accessories		Set includes 2 AAA batteries (LR03) for the monitor, Instruction Manual, Storage case

- * Applied part for this device is Cuff.
- * A range in barometric pressure 700 hPa to 1060 hPa
- * Overall system accuracy of this device is met with the requirement of item 7.9 of EN1060-3.
EN1060-3 : Non-invasive sphygmomanometers Part 3 (European Standard)

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Copyright © Diasense Korea Co.,Ltd. All rights reserved.

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